


产品说明书

CD11c Antibody

货号: ABB5206



产品名称	CD11c Antibody
来源宿主	Rabbit
反应种属	Human
克隆类型	Monoclonal
同种型	IgG
标记	unconjugated
纯化方式	Affinity-chromatography
形式	Liquid
存储溶液	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
存储方式	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
应用	WB IHC-P/mlHC ICC/IF IP
使用方法	WB 1:500~1:2000 IHC-P/mlHC 1:50~1:200 ICC/IF 1:50~1:200 IP 1:50
有效期	one year
别名	p150/95 integrin alpha chain; Leu M5; SLEB6; ITGAX; CD11c; p150; Myeloid membrane antigen alpha subunit; Leukocyte surface antigen p150 95 alpha subunit; Integrin alpha-X;
免疫原	A synthesized peptide derived from human CD11c
SwissProt	P20702
分子量	128kDa
产品介绍	Integrin alpha-X/beta-2 is a receptor for fibrinogen. It recognizes the sequence G-P-R in fibrinogen. It mediates cell-cell interaction during inflammatory responses. It is especially important in monocyte adhesion and chemotaxis.

For research use only!

慧蓝生物 竭诚为您服务

地址: 上海市浦东新区周浦镇天雄路588弄, 电话: 19101712317, 邮箱: 469340997@qq.com