

产品说明书

 α -SMA

货号: ABB00025

产品名称	α -SMA
来源宿主	Rabbit
反应种属	Human Mouse Rat
克隆类型	Rabbit monoclonal
克隆号	S142
同种型	IgG
标记	unconjugated
纯化方式	Protein A affinity purified.
形式	Liquid
存储溶液	PBS (pH7.4), 0.05% BSA, 40% Glycerol. Preservative: 0.05% Sodium Azide.
存储方式	Store at +4°C after thawing. Aliquot store at -20°C. Avoid repeated freeze / thaw cycles.
应用	IHC-P, mIHC, WB, IF
使用方法	WB 1:200 IHC-P 1:250-1: 500 mIHC 1:250-1: 500 IF 1:25
有效期	one year
别名	Alpha-actin-2, ACTA2, Smooth muscle actin, actin-sm, ASM-1, CGA7, ASMH12
SwissProt	P62736
细胞定位	Cytoskeleton, Cytoplasm
分子量	Predicted band size: 42 kDa
产品介绍	Alpha smooth muscle actin, encoded by ACTA2, is expressed in abundance in vascular smooth muscle, comprising 50–70% of total actin, with the remainder composed of β -cytoplasmic and γ -actins. Whereas α -smooth muscle expression is normally restricted to smooth muscle cells, it can also be expressed in non-muscle cells, most notably myofibroblasts that use cell traction forces to remodel extracellular matrix.

For research use only!

慧蓝生物 竭诚为您服务

地址: 上海市浦东新区周浦镇天雄路588弄, 电话: 19101712317, 邮箱: 469340997@qq.com