


# 产品说明书

## CD4 (Y0610) Rabbit mAb

货号: HL18419



产品名称	CD4 (Y0610) Rabbit mAb
来源宿主	Rabbit
反应种属	Human, Mouse
克隆类型	Rabbit monoclonal
克隆号	Y0610
同种型	IgG
标记	unconjugated
纯化方式	Protein A affinity purified.
形式	Liquid
存储溶液	PBS (pH7.4), 50% Glycerol, 0.05% BSA, 0.05% Proclin 300
存储方式	Store at +4°C after thawing. Aliquot store at -20°C. Avoid repeated freeze / thaw cycles.
应用	WB, IHC, mIHC, IF
使用方法	IHC 1:100-1:500; mIHC 1:100-1:500; WB 1:1000-1:5000; IF 1:100-1:500
有效期	one year
别名	T-cell surface glycoprotein CD4; T-cell surface antigen T4/Leu-3; CD antigen CD4;
SwissProt	P01730
细胞定位	Membrane
分子量	60kD
产品介绍	CD4 molecule(CD4) Homo sapiens This gene encodes a membrane glycoprotein of T lymphocytes that interacts with major histocompatibility complex class II antigens and is also a receptor for the human immunodeficiency virus. This gene is expressed not only in T lymphocytes, but also in B cells, macrophages, and granulocytes. It is also expressed in specific regions of the brain. The protein functions to initiate or augment the early phase of T-cell activation, and may function as an important mediator of indirect neuronal damage in infectious and immune-mediated diseases of the central nervous system. Multiple alternatively spliced transcript variants encoding different isoforms have been identified in this gene. [provided by RefSeq, Aug 2010],

For research use only !

**慧蓝生物 竭诚为您服务**

地址: 上海市浦东新区周浦镇天雄路588弄, 电话: 19101712317, 邮箱: 469340997@qq.com