

产品说明书

ZEB2 Rabbit pAb

货号: HLD2733

产品名称	ZEB2 Rabbit pAb
来源宿主	Rabbit
反应种属	Human Mouse
克隆类型	Rabbit polyclonal
同种型	IgG
标记	unconjugated
纯化方式	Protein A affinity purified.
形式	Liquid
存储溶液	PBS (pH7.4),0.05% BSA,40% Glycerol. Preservative: 0.05% Sodium Azide.
存储方式	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.
应用	WB,IHC
使用方法	WB 1:500-1:1000; IHC 1:50-1:100
有效期	one year
别名	SIP-1; SIP1; Smad interacting-protein 1; Smad-interacting 1; SMADIP1; ZEB2; ZFX1B; zinc finger E-box-binding 2; zinc finger homeobox 1b
免疫原	A synthesized peptide derived from human ZEB2
SwissProt	O60315
细胞定位	Nucleus,Chromosome
分子量	Predicted band size: 136 kDa
产品介绍	Transcriptional inhibitor that binds to DNA sequence 5'-CACCT-3' in different promoters. Represses transcription of E-cadherin.Represses expression of MEOX2.

For research use only!

慧蓝生物 竭诚为您服务

地址: 上海市浦东新区周浦镇天雄路588弄, 电话: 19101712317, 邮箱: 469340997@qq.com