

产品说明书

CDX-2 Mouse mAb (Y2)

货号: HL06305

产品名称	CDX-2 Mouse mAb (Y2)
来源宿主	Mouse
反应种属	Human
克隆类型	Mouse monoclonal
克隆号	Y2
同种型	IgG
标记	unconjugated
纯化方式	Protein A affinity purified.
形式	Liquid
存储溶液	PBS (pH7.4), 50% Glycerol, 0.05% BSA, 0.05% Proclin 300
存储方式	Store at +4°C after thawing. Aliquot store at -20°C. Avoid repeated freeze / thaw cycles.
应用	IHC, mlHC, IF
使用方法	IHC 1:100-1:200; mlHC 1:100-1:200; IF 1:50-1:100
有效期	one year
别名	CDX2; CDX3; Homeobox protein CDX-2 (CDX-3) (Caudal-type homeobox protein 2)
SwissProt	P47902; Q99626; O14627
细胞定位	Nuclear
产品介绍	This gene is a member of the caudal-related homeobox transcription factor gene family. The encoded protein is a major regulator of intestine-specific genes involved in cell growth and differentiation. This protein also plays a role in early embryonic development of the intestinal tract. Aberrant expression of this gene is associated with intestinal inflammation and tumorigenesis.

For research use only !

慧蓝生物 竭诚为您服务

地址: 上海市浦东新区周浦镇天雄路588弄, 电话: 19101712317, 邮箱: 469340997@qq.com